

xploris

CODING

My thermometer

xploris

CODING

MY THERMOMETER


- 1 Introduction
- 2 Activity setup
- 3 Coding
- 4 Questions
- 5 Activity summary

1 Introduction

I'm sure you've played a video game or used a smartphone at some point. **Have you ever noticed how they tend to heat up after a while?**


Both people and devices perform best when they're at the right temperature. If it gets too hot, we feel drained and struggle to focus. Devices work the same way—they can slow down, malfunction, or even get damaged if they overheat. **It's like a runner trying to race under the burning sun!**

In this activity, you are going to become a scientist and a programmer. You will use block language and the Xploris device in conjunction with the external temperature probe to build your own thermometer that will let you know if the temperature is comfortable for people.



2 Activity setup

Make sure you connect the external temperature probe to the back of the Xploris as shown in the image.



2

Activity setup




Turn on your Xploris and connect it to your computer or tablet.



Open the XploriLab software on your computer or tablet.




Once inside XploriLab, select the icon to connect the device via USB cable or bluetooth as applicable.



Go to the ENGINEERING section and then to CODING.

3 Coding

In the Coding window, you will find the tools you need to create code using blocks.



3 Coding

The available tools represent blocks that allow you to perform various actions.

Functions

OUTPUT	Output blocks to activate Xploris switches, power source and/or servos.	LOOPS	Loop blocks to perform an action continuously or as long as their conditions are met.
DISPLAY	Blocks to control the Xploris screen.	MATH	Mathematical blocks, such as +,- and many other math functions.
SPEAKER	Blocks to control the Xploris speaker: play sound tracks, notes and control the speaker volume.	VARIABLES	Blocks for creating variables, assign and replace their values.
INPUT	Blocks that enables you to use all Xploris keys and sensors, such as temperature, light, distance, sound and voltages.	LOGIC	Logic operators that will allow decisions to be made based on the state of the data.

3 Coding

1



The screenshot shows the Xploris programming interface. On the left is a vertical menu with categories: OUTPUT (orange), DISPLAY (green), SPEAKER (yellow), INPUT (blue), LOOPS (pink), and MATH (red). The 'DISPLAY' category is selected and highlighted in green. A red checkmark is placed over the 'DISPLAY' label. On the right, a workspace contains several blocks: a 'scroll screen' block (green) with 'up' selected and 'step' set to '10'; a 'draw rectangle' block (green) with '1' and '256' in input fields; a 'display text' block (green) with the text 'Hello Xploris!' and 'position' set to 'center'; and a 'display variable' block (green) with 'position' set to 'center'. A small red and white thermometer icon is visible at the bottom of the workspace.

In this activity, you'll blend the art of drawing with the magic of programming, using the Xploris device and block programming. Your goal is to create a thermometer that displays a specific temperature range. To achieve this, you will begin by exploring the **DISPLAY** function.

2

You will use the **clear screen** block of the **DISPLAY** group to clear the screen of the device, ensuring that you have a clean space to display the results of your programming. To use it, select it and move it to the workspace on the right.

3 Coding

3




Search in **DISPLAY** and drag  under the previous block, remember to join them!

 Set the draws rectangles 3 parameters :

- **First and second parameters** are where oes the rectangle begin and end.
- **Thirdparameter** is color .

4




Let's paint the whole screen blue! You will use the block used to draw rectangles to cover the whole screen, from the top corner to the bottom corner!

Use for the start pixel and **256** for the end pixel. And don't forget the color: **blue**.


3 Coding

5



The Xploris screen is made up of **256 dots**, known as **pixels**. To count them, you start at the top left corner, move to the right and, when you reach the end of the row, move down to the next row and start again from the left.

6



Continue drawing the thermometer. To do this, you will use the colors black, red and white on the Xploris screen. First, drag the `draw rectangle 1 256` option twice, one below the other. In the first rectangle set the **start pixel to “23”** and the **end pixel to “170”**, and choose the color black. In the second rectangle, set the **start pixel to “166”** and the **end pixel to “235”**, also with black color. Thus, you will have drawn two rectangles that will form the outline of our thermometer.

3 Coding

7

```

clear screen
draw rectangle 1 256 blue
draw rectangle 23 170 black
draw rectangle 166 235 black
draw rectangle 40 169 white ✓
draw rectangle 183 218 red ✓
    
```

To draw the temperature indicator in red and white, drag two more `draw rectangle 1 256 blue` and join it to the ones previously added. In the first one, set the start point to “40” and the end point to “169”, and choose the color white. In the second one, set the start point to “183” and the end point to “218”, and choose the color red. This way you have drawn the base of your thermometer

8

In the **LOOPS** group section, use the `forever` block.

This will let you repeat the instructions inside it endlessly. **Simply place it next to the work you've already completed!**

Now, we are going to program the feature that lets us see how the temperature changes. We'll be working within a range of **20 to 30 degrees Celsius.**

3 Coding

9

- From **VARIABLES**, you are going to use the **set y to** block to create a “y” variable that will have the value of the temperature, which will be between a comfortable range for it to work well.
- Add it inside the **forever** block.
- Then, in order to operate the variable more comfortably, from group **MATH** drag **+** inside “to”, to the interior of **set y to**.

```

clear screen
draw rectangle 1 256
draw rectangle 23 170
draw rectangle 166 235
draw rectangle 40 169
draw rectangle 183 218
forever
  set y to +
  
```

3 Coding

10

- Now you are going to fill in the information requested by the operator block .
- Go to group **INPUT** and drag the block **temperature celsius** and from **MATH** drag the **0** numeric block.
- Then fill the left side of the  operator block with **temperature celsius** and make sure the scale you choose is **“Celsius”**.
- On the right side of the operator block, add the numeric block and change the value to **“20”**. The  operator should be configured to act as a **“-” (subtraction)**.

Based on what you did before, you will be giving the variable “y” the value that the Xploris device reads from the **external temperature probe**, together with the instruction to subtract 20 from this value .

```

clear screen
draw rectangle 1 256
draw rectangle 23 170
draw rectangle 166 235
draw rectangle 40 169
draw rectangle 183 218
forever
  set y to temperature celsius - 20
  
```


3 Coding

11

Continuing with the programming:

- Go to the group **LOGIC** and drag the function  just below what you have already done.
- Then you will use 3 more functions from the same group: one , and you add it inside the previous function, and two  and you drag them to the workspace.
- From the group **VARIABLES** you will take two  and from the group **MATH** you will use two numeric blocks .

When you take the programming blocks to place them inside other blocks, make sure you leave them right up against the left edge of the block where you want to place them. If you don't do it like this, they may not fit inside the block correctly.




3


Coding

12

Now it's time to set the temperature at which your thermometer is going to make changes. To do this, you will have to indicate the values between which you want it to work.

- Make sure that both variables  have the letter 'y' assigned to them and enter the value '0' in the numeric block  on the left and the value '9' in the one on the right.
- Then, in the left comparator  use the symbol " \geq " and in the right comparator the symbol " $<$ ".
- Don't forget to set the block  to "and"!

That's it, you're almost done!



```


forever
  set y to temperature celsius - 20
  if y >= 0 and y < 9
    do
  
```

3 Coding

13

Now, you'll add color to your thermometer by following these steps, but only if the previously specified conditions are met.

- From the group **DISPLAY**, drag two **draw rectangle** blocks and put them inside the **if** block in the "do" part.
- Then, from **MATH**, drag two **+** blocks and two **0** numeric blocks.
- Finally, from **VARIABLES**, drag the **i** variable and use it in the workspace.



3 Coding

14

You're almost done! Now it's time to put the information inside the blocks.

- For the first `draw rectangle` block, the start position will be at pixel "40" and the end position at pixel "169", both intend white.
- In the second `draw rectangle` block, already assembled, fill in the information from left to right: put a numerical block with the value "168", then the `+` block with the subtraction option, and in the block on the right, put the variable "selected in "y" doing `i`, the operation "x" (**multiplication**) assigning it the value "16".
- Then fill the `draw rectangle` block with the final pixel '169' and make it **red**.

```

forever
  set y to temperature celsius - 20
  if y >= 0 and y < 9
    do
      draw rectangle 40 169 white ✓
      draw rectangle 168 - y x 16 169 red ✓
  
```

3 Coding

14

- Finally, from the **LOOPS** group use the **delay 100 (ms)** block, which will help you to see the temperature variation on the screen of the Xploris device more clearly.

```

clear screen
draw rectangle 1 256
draw rectangle 23 170
draw rectangle 166 235
draw rectangle 40 169
draw rectangle 183 218
forever
  set y to temperature celsius 20
  if y >= 0 and y < 9
  do
    draw rectangle 40 169
    draw rectangle 168 y * 16 169
  delay 100 (ms)
  
```


3 Coding

To make sure that the program works correctly, we will follow these final steps:

Press the three-bar icon at the top and select the “Save” option. Then, assign a name and save our program. 

Press the “Upload” button in the Xplorilab interface. This will transfer the program to the Xploris device. 

Once the program is loaded, press the “Play” button on Xplorilab software. Take a moment to admire your work: an amazing thermometer brought to life through block programming. 



The screenshot shows the Xploris planet interface with a cloud icon and the text "Xploris planet". Below this, there are two buttons: "Upload" and "Open". A horizontal line separates this section from the next. The next section has a laptop icon and the text "Local". Below this, there are two buttons: "Save" (with a red checkmark) and "Open". A final horizontal line separates this from the last section, which has a book icon and the text "Lesson Plans".



4

Questions

1

Sciences

How could we use this code to represent other scientific data, such as distance or the amount of light?

2

Art

What patterns or shapes could we create if we changed the values of the rectangles?

3

Let's keep experimenting!

Could we get the program to emit sounds as the temperature varies?

5

Activity summary



We used the Xploris screen to display a thermometer that tracks temperature variations



We used the values from a temperature probe and saved them in a variable. We used a “Logic” programming block so that the program only draws the elements when the temperature is within a certain range.



We used the XploriLab program to draw a thermometer on the screen that adjusts based on the temperature measured by one of the Xploris device probes.



xploris
CODING

My thermometer