



xploris

SCIENCE

How loud: I can't hear you!

xploris

SCIENCES

I CAN'T HEAR YOU!

- 1 Introduction
- 2 Setting up the experiment
- 3 Data collection
- 4 Data analysis
- 5 Questions
- 6 Activity summary

1 Introduction

Have you ever seen two people talking to each other and wanted to hear what they were saying, but couldn't? Either because they are talking too quietly or because they are too far away, sometimes we simply can't hear the sound.

In this lesson you will study the relationship between sound and distance, using the Xploris microphone sensor.

The question you will answer will be:



How does the sound change as I get closer to the source?

2 Setting up the experiment



Turn on your Xploris and connect it to your computer or tablet.



Open the XploriLab software on your computer or tablet.



Once inside XploriLab, select the icon to connect the device via cable or bluetooth as applicable.




Go to the SCIENCE section and then to DATA LOGGER.

2

Setting up the experiment

Find a quiet hallway in your school. Place a speaker at one end and move away from it by marking distances on the floor separated by one meter until you reach 5 meters. Then, one person should turn on the speaker at high volume and you should register the sound with the Xploris sensor as you get closer to the speaker.



It is necessary for the speaker to produce a clear single frequency sound.

To do this, use any online tone generator and select a frequency close to 1000 Hz.



2

Setting up the experiment

↖ XploriLab software configuration

1

To start any configuration related to the sensors, please select the “setup” icon.

The sensor you will use for this activity is the **sound** sensor and you will configure it to take **manual measurements with a total of 10 samples**.

Once the configuration has been completed, select “Apply” to save it.



Choose Sensor ✕

Light

Voltage left Voltage Right

Ext. Temperature Amb. Temperature

Distance Speed

Pulse Heart rate

Sound

Rate Samples

Manual 10

Apply




2

Setting up the experiment

↖ XploriLab software configuration

1

Click the “line” icon on display and select the “bars” option.




2

Click the PLAY icon to start recording. Then, when you want to collect a new measurement - press the Sound key on the Xploris.



3 Data collection

Once you have made the marks on the floor every one meter, turn on the speaker as loud as possible and measure the sound manually as you get closer to the speaker. Remember to make sure that the place is as quiet as possible so as not to interfere with the measurements.




4 Data analysis

1

Use markers to add labels to the columns of the graph. To do this you must select the "Marker" icon:



GRAPH WITH MARKERS




4

Data analysis

2

You can add photos to the notes within a graph by doing the following:

1. Select the note icon. 
2. Click on the bar where you want to add a note.
3. A dialog box will open and allow you to add the note with text and images.


Write your note



4 Data analysis

3

GRAPH WITH NOTES



5

Questions

1

Let's take a look at the graph

What was the highest sound intensity detected and what was the lowest?

2

Let's evaluate the data

What relationship can be made between sound intensity and distance from the speaker?

3

Let's investigate!

Look up in books or on the internet how we can hear and what happens to the ears when we are exposed to very loud sounds.

4

Let's keep experimenting!

What do you think would happen to the measurements if instead of having the speaker at moderate volume we use it at very high volume? Formulate your hypothesis and test it using the Xploris microphone sensor.

6

Activity summary



We use the Xploris sound sensor to study the sound intensity produced by a speaker.



We analyzed the data to establish the relationship between distance from the speaker and sound intensity.



We answered questions by analyzing the values obtained, as well as investigating how hearing occurs and what happens when we are exposed to loud sounds. Finally, we changed the volume of the speaker and repeated the experiment to see if there were any differences.



xploris
SCIENCES

How loud: I can't hear you!