

Applied Sciences

▶ Concentrations that absorb

Measuring percentage transmittance of solutions at different concentrations



Digital content
provided by
 **efecto educativo**
efectoeducativo.com

FRAMEWORK FOR K-12 SCIENCE EDUCATION © 2012

The Dimension I practices listed below are called out as **bold** words throughout the activity.

Dimension 1 Science and Engineering Practices	✓	Asking questions (for science) and defining problems (for engineering)	✓	Use mathematics and computational thinking
	✓	Developing and using models	✓	Constructing explanations (for science) and designing solutions (for engineering)
	✓	Planning and carrying out investigations	✓	Engaging in argument from evidence
	✓	Analyzing and interpreting data	✓	Obtaining, evaluating, and communicating information
Dimension 2 Cross Cutting Concepts	✓	Patterns	✓	Energy and matter: Flows, cycles, and conservation
	✓	Cause and effect: Mechanism and explanation		Structure and function
	✓	Scale, proportion, and quantity		Stability and change
	✓	Systems and system models		

Applied Sciences

Concentrations that absorb

Measuring percentage transmittance of solutions at different concentrations

USA Standards Correlation

Dimension 3 Core Concepts	Discipline	Core Idea Focus
	Physical Science	PS1: Matter and Its Interactions
		PS1.A: Structure and Properties of Matter

NGSS Standards	Middle School Standards Covered	High School Standards Covered
	MS.PS-SPM: Structure and Properties of Matter	HS.PS-SPM: Structure and Properties of Matter

NATIONAL SCIENCE EDUCATION STANDARDS © 2002

Content Standards (K-12)			
✓	Systems, order, and organization		Evolution and equilibrium
✓	Evidence, models, and explanation	✓	Form and Function
✓	Constancy, change, and measurement		

Physical Science Standards Middle School		Physical Science Standards High School	
✓	Properties and Changes of Properties in Matter		Structure of Atoms
	Motions and Forces	✓	Structure and Properties of Matter
	Transfer of Energy		Chemical Reactions
			Motions and Forces
			Conservation of Energy and Increase in Disorder
			Interactions of Energy and Matter

LEARNING OBJECTIVES

Core Objectives (National Standards):

- Develop the ability to refine ill-defined questions and direct to phenomena that can be described, explained, or predicted through scientific means.
- Develop the ability to observe, measure accurately, identify and control variables.
- Decide what evidence can be used to support or refute a hypothesis.
- Gather, store, retrieve, and analyze data.
- Become confident at communicating methods, instructions, observations, and results with others.

Activity Objectives:

The purpose of this activity is to relate light absorbance and transmittance in solutions with different concentrations, create a hypothesis and proceed to test it using the Globisens Labdisc colorimeter sensor.

Time Requirement:

60 - 90 Minutes

Applied Sciences

Concentrations that absorb

Measuring percentage transmittance of solutions at different concentrations

Objective

The purpose of this activity is to relate light absorbance and transmittance in solutions with different concentrations, create a hypothesis and proceed to test it using the Labdisc colorimeter sensor.

Applied Sciences

Concentrations that absorb

Measuring percentage transmittance of solutions at different concentrations

Introduction and theory

The aim of the introduction is to focus students on the lesson subject by refreshing acquired knowledge and asking questions which encourage research development. Key concepts from the theoretical framework, applied by the students during the lesson, are taught.

Introduction

As part of our everyday life we prepare and use different solutions. Each one is composed of two parts: The solvent (most of the time water), i.e. the substance into which the solute is dissolved; and the solute, which is the substance that is dissolved in the solvent. A common example of a solution may be salt water for cooking pasta, which becomes more salty the more concentrated the solution is.

When you are doing scientific work, particularly in chemistry, biology or medical sciences, it is very important to know exactly the concentration of the solutions you are working with. To measure concentration, several techniques can be used.

Applied Sciences

Concentrations that absorb

Measuring percentage transmittance of solutions at different concentrations

Introduction and theory

?

Can you think of some different solutions you have prepared at school or at home?

?

When you are making tea, how do you know when it is ready? Describe.

Carry out the experiment activity with your class so that at the end you'll be able to answer the following question:

?

How can a light beam that passes through a sample of a given solution help determine its concentration?


Theoretical

A solution is a homogeneous mixture of two or more substances that exist in a single concentration. To determine the concentration of a given solution you can apply the **Beer-Lambert law**. This law states that the concentration of a solute is proportional to the absorbance.

The colorimeter allows light to pass through a cuvette containing a sample of the solution which absorbs some of the incoming beam. When the ray of light of a given wavelength and intensity I comes perpendicularly into contact with the solution of a tinted chemical compound, the compound will absorb part of the radiation (I_a). The remaining part of the solution (I_b) will strike a detector. As such the following equation is demonstrated:

$$I_0 = I_a + I_b$$

The absorption of light is related to the number of molecules present in the solution (concentration of the solution).



The Beer-Lambert law defines the relationship between the concentration of a solution and the amount of light absorbed by the solution:

$$A = \epsilon dC$$

Where:

A = Absorbance

ϵ = Molar absorptivity [$\text{L mol}^{-1} \text{cm}^{-1}$]

d = Path length of the cuvette containing the sample [cm]

C = Concentration of the compound in solution [mol L^{-1}]

Transmittance is the relationship between the amount of light that is transmitted to the detector once it has passed through the sample (I) and the original amount of light (I_0). This is expressed in following formula:

$$T = \frac{I}{I_0}$$

Where I_0 is the intensity of the incident light beam and I is the intensity of the light coming out of the sample. Transmittance is the relative percent of light passed through the sample. Thus, if half the light is transmitted, we can say that the solution has a 50% transmittance.

$$T\% = \frac{I}{I_0} \cdot 100\%$$

The relations between **transmittance** (T) and **absorbance** (A) can be expressed by following:

$$A = \log_{10}\left(\frac{1}{T}\right)$$

$$A = \log_{10}\left(\frac{100}{T[\%]}\right) = 2 - \log T[\%]$$

Applied Sciences

Concentrations that absorb

Measuring percentage transmittance of solutions at different concentrations

Introduction and theory

Now students are encouraged to raise a hypothesis which must be tested with an experiment.

?

If you have one solution of different concentrations, how do you expect the transmittance percentages to change when you measure from the lowest to the highest concentration? Why?

Applied Sciences

Concentrations that absorb

Measuring percentage transmittance of solutions at different concentrations

Activity description

Students will study the relationship between transmittance, absorbance and concentration of one solution using the Lambert-Beer law. They will calculate the concentration of a sample using mathematical tools for graph analysis.

Preparation of the samples

The teacher will prepare two samples of solution for the students. The students will obtain the concentrations of the samples through the experiment.

The instructions to prepare the samples are as follows:

	SAMPLE 1	SAMPLE 2
Solute (instant coffee) [g]	1	1
Solvent (distilled water) [mL]	140	90
Concentration [g/mL]	0.0071	0.011

The teacher will give 3 mL. of one of the samples to each work group. The students will be told that they have to obtain the concentration. At the end of the class they will compare their results with the theoretical value.

Applied Sciences

Concentrations that absorb

Measuring percentage transmittance of solutions at different concentrations

Activity description

Students will read a graph to obtain the concentration. To understand it they have to know the slope-interception form for the equation of a line.

$$Y = mX + n$$

Where

m = slope

X = x coordinate



Y = y coordinate

n = y intercept (where the line crosses the y axis)

Students must:

1. Obtain the transmittance percentage by using the colorimeter.
2. Divide the transmittance percentage by 100 to obtain transmittance (T) of the sample.
3. Calculate absorbance (A) using the following formula: $A = \text{LOG}_{10}T$.
4. Obtain X by solving the slope-interception equation, substituting the y coordinate for the activity description absorbance value.

- 1 Labdisc
- 2 USB connector cable
- 3 Cuvette for the colorimeter
- Distilled water
- 100 mL beaker
- Instant coffee
- Paper towel
- 4 test tubes
- Wash bottle
- Stirrer bar
- Balance



Applied Sciences

Concentrations that absorb





Measuring percentage transmittance of solutions at different concentrations


Resources and materials




a. Labdisc Configuration

To collect measurements with the microphone sensor, the Labdisc must be set up according to the following steps:

- 1 Open the GlobiLab software and turn on the Labdisc.
- 2 Click on the Bluetooth icon in the bottom right corner of the GlobiLab screen. Select the Labdisc you are using currently. Once the Labdisc has been recognized by the software, the icon will change from a grey to blue color   2/127 . If you prefer a USB connection follow the previous instruction clicking on the USB icon. You will see the same color change when the Labdisc is recognized   0/127 .

- 3 Click on  to configure the Labdisc. Select colorimeter in the “Logger Setup” window. Enter “Manual” for the sample frequency and the number of samples.





Applied Sciences

Concentrations that absorb

Measuring percentage transmittance of solutions at different concentrations

Using the Labdisc

- 4 Once you have finished the sensor configuration start measuring by clicking .
- 5 Once you have finished measuring stop the Labdisc by clicking .

Applied Sciences

Concentrations that absorb

Measuring percentage transmittance of solutions at different concentrations

Experiment

The following steps explain how to perform the experiment:

Note: Before you start recording, calibrate the colorimeter with the solvent you will use during the experiment. This step is important to minimize the experimental error of the measurement. To calibrate, add some solvent (in this case, distilled water), filling $\frac{3}{4}$ of the sensor's cuvette and place the cuvette in the output. Then, press the button on the colorimeter until you hear the Labdisc's signal. Once you have finished the calibration remove the cuvette from the colorimeter.

- 1 Mark each beaker with the number of the sample (1 to 6).
- 2 Place distilled water in beaker 1.

- 3 Prepare samples 2 to 6 as shown in the table below by measuring the amount of distilled water (solvent) and adding 1 gram of coffee (solute) to each glass.

	Sample 2	Sample 3	Sample 4	Sample 5	Sample 6
Coffee [gr]	1	1	1	1	1
Distilled water [mL]	150	120	100	80	60
Concentration [g/mL]	0,0066	0,0083	0,0100	0,0125	0,0166

- 4 Use the stirrer bar to mix the coffee until it dissolves completely.
- 5 Measure the transmittance of the solution. To take a data reading, fill the cuvette until it is $\frac{3}{4}$ full with solution from sample 1. ALWAYS hold the upper part of the cuvette to ensure the sides remain clean from fingerprints.
- 6 Clean and dry the cuvette with paper towels.


Applied Sciences

Concentrations that absorb

Measuring percentage transmittance of solutions at different concentrations

Experiment

- 7 Place the cuvette in the output of the colorimeter. Cover the cuvette moving the edge of the Labdisc and start measuring.



Applied Sciences

Concentrations that absorb


Measuring percentage transmittance of solutions at different concentrations

Experiment

- 8 Record the transmittance percentage. Take the sample out of the sensor and clean the cuvette with distilled water. Repeat the steps with the other samples.
- 9 Once you have finished measuring stop the Labdisc.

The graph below should be similar to the one the students came up with. The decreasing order of percentages can also be observed in both green and blue colorimeters.

Transmittance percentage as a function of concentration:




Applied Sciences

Concentrations that absorb

Measuring percentage transmittance of solutions at different concentrations

Results and analysis


The following steps explain how to analyze the experiment results:

- 1 Export the data to Excel clicking on  . Save the data on your computer.
- 2 Observe the percentage transmittance values obtained when we illuminate the samples with different light beam colors (different wavelengths: Red, green and blue).

- 3 Organize the data and calculate transmittance and absorbance values starting from the transmittance percentage given by the Labdisc. Data should be organized as shown in the following table:

Sample – Concentration [gr/mL] - %Transmittance red – Transmittance red – Absorbance red - %Transmittance green – Transmittance green – Absorbance green - %Transmittance blue – Transmittance blue – Absorbance blue

	A	B	C	D	E	F	G	H	I	J	K
1	Muestra	Concentración [g/mL]	%Transm. rojo	Transm. rojo	Abs. rojo	%Transm. Verde	Transm. Verde	Abs. Verde	%Transm. Azul	Transm. Azul	Abs. Azul
2	1	0		=C2/100	=-LOG10(D2)		=F2/100	=-LOG10(G2)		=I2/100	=-LOG10(J2)
3	2	0,0066									
4	3	0,0083									
5	4	0,01									
6	5	0,0125									
7	6	0,0166									

- 4 Label column C as percentage transmittance of red, column F as percentage transmittance of green and column I as percentage transmittance of blue. Use the formula of the table to calculate the transmittance (columns D, G and J) and absorbance (columns E, H and K) of each sample.
- 5 Observe and compare the transmittance values for each sample. Remember that the Lambert-Beer law can be used for absorbance values between 0 and 1. Therefore, select the color where you obtained most results within this absorbance range.
- 6 Create a line graph of absorbance as a function of concentration. To do this, first create a scatter plot by clicking on the  **XY (Scatter)** button.
- 7 Select the Data Range. Select the concentrations value for the x axis and absorbance values for the y axis.
- 8 Add a regressing line to the graph by right-clicking on the points and selecting “Add Trendline”. Select the linear regression type. Choose the options chart and select “Display equations on chart”. Press Ok to finish.

Applied Sciences

Concentrations that absorb

Measuring percentage transmittance of solutions at different concentrations

Results and analysis

?

How do the results relate to your initial hypothesis? Explain.

?

What was the relationship between the transmittance and concentration of the solutions?

?


What was the relationship between the absorbance and concentration of the solutions?

The graph below should be similar to the one the students came up with.

Concentration as a function of absorbance:

Y axis: Absorbance [red colorimeter]

X axis: Concentration [gr/mL]



Applied Sciences

Concentrations that absorb

Measuring percentage transmittance of solutions at different concentrations

Conclusions

Following are some questions and answers which should be developed by students in order to elaborate on their conclusions.

?

What variables are related in the Lambert-Beer law?

Students should point out that the Lambert-Beer law relates absorbance of a solution to its concentration.

?

If you have a solution and samples of it with different concentrations; in which range of absorbance can we use the Lambert-Beer law?

Students should go over the information given in the theoretical framework and point out that the Lambert-Beer law is only valid for absorbance values between 0 and 1.

Applied Sciences

Concentrations that absorb

Measuring percentage transmittance of solutions at different concentrations

Conclusions

?

Suppose the equation on the Excel graph of absorbance as a function of concentration is a straight line. What is the relationship between the expression and the variables related in the Lambert-Beer law?

Students should understand the meaning of the equation of the graph. In this case, the y coordinate represent absorbance and the x coordinate concentration. Therefore, we can understand the equation as following expression:

$$\text{Absorbance} = m * \text{Concentration} + n$$

?

If you have a sample of unknown concentration, how would the Lambert-Beer law be useful?

Students should point out that you can obtain the concentration of an unknown sample by measuring absorbance and using the Lambert-Beer's law to compare its absorbance to a known solution as an absorbance reference.

Applied Sciences

Concentrations that absorb

Measuring percentage transmittance of solutions at different concentrations

Activities for further application

The aim of this section is for students to extrapolate the acquired knowledge during this class through its application in different contexts and situations. Furthermore, it is intended that students question and present possible explanations to the experimentally observed phenomena.

Further questions:

?

Suppose you have an instant coffee solution of an unknown concentration which you want to calculate using the Lambert-Beer law. However, your absorbance is higher than one - how can you calculate concentration?

Students should point out that they will have to decrease the absorbance value to use the Lambert-Beer law, diluting the initial sample to obtain an absorbance value between 0 and 1. Starting from that point, they could calculate the concentration of the diluted solution and then the concentration of the original sample.

Applied Sciences

Concentrations that absorb

Measuring percentage transmittance of solutions at different concentrations

Activities for further application

?

What was the concentration of the sample prepared by your teacher?

Students should use the knowledge obtained during the class to calculate the concentration of the work sample and obtain 0.0071 [gr/mL] or 0.011 [gr/mL] approximately, depending on the sample they worked with.

efectoeducativo

globisens



 **Globisens**
Lab classes with sensors
Labdisc

 Digital content
provided by
 efecto educativo
efectoeducativo.com